

# bliss fundamental commands

## Starting a BLISS session

Pandora's box

In a terminal: `. blissenv [-d]` (note the dot and space)

`bliss -s <session_name>`

`bliss -h`: help about bliss command

To **detach** from a session: `ctrl-b d`

To **exit** from a BLISS session: `exit` `↵` or `ctrl-d`

## Counters

Counters and Measurement Groups (MG)

`lscnt()`: print list of all counters

`lsmg()`: print list of all measurement groups

`ACTIVE_MG`: info about the current measurement group  
Measurement groups are created in config file.

`menu(myMG)`: open dialog to edit **myMG**

`myMG.set_active()`: define **myMG** as default M.G.

`ct(0.2, i0diode)`: count for 0.2 second using  
i0diode counter

`sct(0.2, i0diode)`: `ct` saved as a scan

## Axis

`wa()`: display positions (user & dial) of all motors

`wm(mot1,...motN)`: display positions and limits  
for specified motors

`mv/move(m1, 2.0)`: move motor **m1** to 2.0

`umv(mot1,8)`: move + display position during move

`umvr(mot1, 0.1)`: idem with a move of 0.1 unit  
relative to current position

`m1.position=3.4`: set user position (⇒chg. offset)

`m1.velocity=1.5`: set velocity (uu/s)

`m1.acceleration=5.0`: set acceleration (uu/s<sup>2</sup>)

`m1.offset`: ≠ between *dial* and *position*

`m1.backlash`

`m1.dial=10.0`: set *dial* position (⇒change user pos.)

`m1.limits=(-10,10)`: set low/high limit to -10, 10

`m1.sign`: direction of the movement relative to dial

`m1.tolerance` `m1.steps_per_unit`

`m1.unit` `m1.encoder` `m1.state`

`position=(sign*dial)+offset`

`uu`: user units (mm, μm, deg...)

*properties in italic are Read-Only*

## Standard scans

Common step by step scans

`ascan(motor, start, stop, intervals,  
count_time, *cnt_args)`

Perform a step\_scan from <start> to <stop>  
counting <intervals> times <count\_time> seconds  
using the given counters or Measurement groups.

`a2scan()..a5scan()`: scan with 2..5 movements

`anscan()`: scan with an arbitrary number of motors

`dscan()..d5scan()`: scan relative to the current pos.

`dnscan()`: relative scan with arbitrary nb. of motors

`amesh(mot1, start1, stop1, interv1,  
mot2, start2, stop2, interv2,  
count_time, *counters)`

Absolute 2D scan on a regular grid

`dmesh()`: relative 2D scan

`lineup(...)`: same as `dscan` then goes to max

`timescan(ctime, *cnt_args)`: endless counts

`loopscan(npoints, ctime, *cnt_args)`  
count <npoints> times

`pointscan(mot, pos_list, ctime)`  
scan over a positions list

`lookupscan([(m1, <pos_list1>)...], ctime)`  
scan over a variable number of motors and positions.

## SCAN\_SAVING

Data Saving and DATA POLICY

`SCAN_SAVING` `↵`: display saving parameters

`newproposal("mr1234")`

`newsample("kryptonite")`

`newdataset("Zn_inclusion")`

## Shutters

Safety shutters and front-end

`bsh1.close()`: close shutter **bsh1**

`fe.open()`: open front-end

`fe.mode="AUTOMATIC" / "MANUAL"`

set front-end in automatic or manual opening mode

## Help

Message in a bottle

`help(command)`: print help about <command>

`last_error()`: details about last error occurred

`last_error(-2)`: details about previous error occurred

`prdef(function)`: print the code and location  
of a function.

BLISS objects have a "in-shell info" feature:

typing its name + `↵` print details about it.

`lsobj("*diod*")`: list all session's objects  
with "diod" in their name

## Shell functions

**F3**: enter *history mode*

In history mode: `space` to select,

`ctrl-o`: once to validate, second time to execute

History is saved per user

**F2**: ptypython configuration (colors etc.)

**F4**: switch shortcut mode

**F5**: switch to/from scan view

**F6**: paste mode

**F7**: typing helper (de)activation

**F8**: set logbook filling from shell **on/off**

`ctrl-r xx` `↵`: search for commands starting  
by "xx" in history

`_XX`: reference to the shell output number XX

## Plotting

display and play with data

`flint()`: launch flint plotting tool

`plotselect(diode1)`: select diode1 for plotting

`SCAN_DISPLAY`: object to configure plotting

`SCAN_DISPLAY.auto=False`: disable plotting

`cen()/goto_cen()`: display middle of *fwhm* (and *fwhm*)  
/ move scanned motor to this position

`peak()/goto_peak()`: idem for max

`com()/goto_com()`: idem for center of mass

*For your favorite commands*